

# In Situ Hybridization

Modified from C. Henry, M. Halpern and Thisse labs  
April 17, 2013

## Table of Contents

<b>Reagents .....</b>	<b>2</b>
<b>AP Buffer.....</b>	<b>2</b>
<b>Developing Solution .....</b>	<b>2</b>
<b>Hybridization buffer .....</b>	<b>2</b>
<b>PBT .....</b>	<b>2</b>
<b>PI Buffer (Blocking solution) .....</b>	<b>2</b>
<b>SSCTween .....</b>	<b>2</b>
<b>Stop Solution.....</b>	<b>3</b>
<b>Riboprobe template preparation (from plasmid).....</b>	<b>4</b>
<b>Riboprobe template preparation (from genomic DNA).....</b>	<b>5</b>
<b>Riboprobe synthesis.....</b>	<b>5</b>
<b>ISH Protocol.....</b>	<b>7</b>
<b>Day 1 .....</b>	<b>7</b>
Embryo preparation and fixation .....	7
Hybridization .....	7
<b>Day 2 .....</b>	<b>8</b>
Washing out the probe .....	8
Antibody binding .....	8
<b>Day 3 .....</b>	<b>9</b>
Develop.....	9
Imaging embryos.....	9

## Reagents

### AP Buffer

Component	Location	Amount (per 50 ml)	Amount (per 250 ml)	Final
1M Tris pH 9.5	RT	5 ml	25 ml	100 mM
0.5 M MgCl <sub>2</sub>	RT	5 ml	25 ml	50 mM
5M NaCl	RT	1 ml	5 ml	100 mM
20% Tween-20	RT	250 µl	1.25 ml	0.1%

### Developing Solution

Component	Location	Amount (per 1 ml)	Amount (per 5 ml)
AP Buffer	RT	1 ml	5 ml
100 mg/ml NBT	Freezer ISH components	4.5 µl	22.5 µl
50 mg/ml BCIP	Freezer ISH components	3.5 µl	17.5 µl

### Hybridization buffer

To a 50 ml conical tube add:

Component	Location	Amount (per 50 ml)	Final
Formamide	4 °C	25 ml	50%
20X SSC	RT	12.5 ml	5X
Heparin (50 mg/ml)	Freezer ISH components	50 µl	50 µg/ml
Yeast tRNA (50 mg/ml)	Freezer ISH components	500 µl	500 µg/ml
20% Tween	RT	250 µl	0.1%
H <sub>2</sub> O	Sterile water	10.8 ml	
1M Citric acid, pH 6.0	RT	460 µl	

Label and date

Store up to 2 months at -20 °C

### PBT

50 ml 1X PBS + 250 µl 20% Tween-20

### PI Buffer (Blocking solution)

Component	Location	Amount (per 10 ml)	Amount (per 50 ml)
PBT	RT	9.6 ml	48 ml
BSA (100 mg/ml)	Freezer ISH components	200 µl	1 ml
Sheep serum	Freezer ISH components	200 µl	1 ml

### SSCTween

50 ml 2X SSC + 250 µl 20% Tween-20 (prevents sticking)

## Stop Solution

<b>Component</b>	<b>Location</b>	<b>Amount (per 10 ml)</b>	<b>Amount (per 50 ml)</b>
1x PBS (pH 5.5)	RT	9.9 ml	48 ml
EDTA (0.5 M)	RT	20 $\mu$ l	100 $\mu$ l
20% Tween20	RT	50 $\mu$ l	250 $\mu$ l

## Riboprobe template preparation (from plasmid)

1. Set up digests of 10 µg of plasmid in 100 µl:

<b>Component</b>	<b>Volume</b>
DNA	Vol. to 10 µg
H <sub>2</sub> O	to final volume of 100 µl (including Buffer and enzyme)
Buffer	10 µl
Enzyme (10u/ µl)	2 µl

2. Incubate at 37°C for 1-2 hours or O/N.
3. Run 1 µl on a gel to check for complete linearization.
4. Phenol:chloroform extract:
  - a. Bring volume up to 200 µl w/H<sub>2</sub>O.
  - b. Add 100 µl Tris saturated phenol.
  - c. Add 100 µl chloroform.
  - d. Vortex.
  - e. Spin 2 minutes at maximum speed at RT.
  - f. Carefully transfer top aqueous layer to a clean tube.
5. Chloroform extract:
  - a. Add 200 µl chloroform.
  - b. Vortex
  - c. Spin 2 minutes at maximum speed at RT.
  - d. Carefully transfer top aqueous layer to a clean tube.
6. Ethanol precipitate:
  - a. Add 1/10 volume 3 M NaOAc (pH 5.2).
  - b. Add 2 volumes 100% EtOH.
  - c. Vortex.
  - d. Spin 15 minutes at maximum speed at 4°C.
  - e. Carefully pour off supernatant.
7. Wash pellet by adding 1 ml 70% EtOH.
  - a. Spin 5 minutes at maximum speed at 4°C.
  - b. Carefully pour off supernatant and air dry, inverted
8. Resuspend in 40 µl TE.
9. Proceed to Riboprobe synthesis or store at -20°C.

## Riboprobe template preparation (from genomic DNA)

### 1. PCR reaction to generate probe synthesis template

Assemble the following in a PCR tube:

- a. Genomic DNA (~200 ng)
- b. ddH<sub>2</sub>O to 25  $\mu$ l final volume
- c. 2.5  $\mu$ l PCR buffer (10x stock)
- d. 0.5  $\mu$ l dNTPs (10 mM)
- e. 1  $\mu$ l 10  $\mu$ M stock forward primer
- f. 1  $\mu$ l 10  $\mu$ M stock reverse primer
- g. 0.5  $\mu$ l Taq Polymerase

Perform a standard PCR reaction (adjust annealing temp if needed):

95°C for 3 min with repetition of  
95°C for 30 sec,  
55°C for 30 sec,  
72°C for 30 sec for 30 cycles, then  
72°C for 10 min.

*Check an aliquot (1  $\mu$ l) on an agarose gel to make sure the PCR product is the expected size.*

## Riboprobe synthesis

1. DIG-labeled probe synthesis (this is for T3, which uses a 5x buffer and requires a DTT supplement) – adjust accordingly if using a different polymerase with different reaction requirements (e.g., some brands include the DTT in a 10x buffer):

PCR reaction (not purified) or prepared template	2 $\mu$ l
DEPC H <sub>2</sub> O	5 $\mu$ l
10x NTP mixture (10mM, ATP, CTP, GTP and 6.5 mM UTP)	2 $\mu$ l
10x DIG-UTP (3.5 mM)	2 $\mu$ l
5xTranscription Buffer	4 $\mu$ l
DTT	2 $\mu$ l
RNase Inhibitor (40 U/ $\mu$ l)	1 $\mu$ l
T3 Polymerase (17 U/ $\mu$ l)	2 $\mu$ l
	(total volume: 20 $\mu$ l)

2. Incubate at 37°C for 2 hours.
3. Remove 0.5  $\mu$ l to confirm by gel electrophoresis (should see RNA *and* DNA template).

4. Add 2  $\mu$ l RNase-free DNase (10 U/  $\mu$ l) to reaction tube and incubate at 37°C for 15 minutes.
5. During DNase incubation, prepare column (GE G-50 column).
6. Break off bottom.
7. Spin at 3K rpm for 1 min to remove liquid.
8. After the DNase incubation, transfer the probe mix over prepared column.
9. Spin for 2 min at 3K rpm **into a fresh tube**.
10. Remove 0.5  $\mu$ l and run on a gel with the pre-DNase treatment product (step 3)
11. Test probe and if it's good, make aliquots and store at -20°C to avoid repeated freezing and thawing.

---

# ISH Protocol

## Day 1

### Embryo preparation and fixation

1. Put up to 10 embryos each in microfuge tubes.
2. Under the hood, add 4% PFA to the tubes – WEAR GLOVES.
3. Incubate for 4h at RT or O/N @ 4°C.
4. **If embryos have pigmentation bleach them before dehydration.**
  - a. Incubate the embryos at room temperature in a 3% H<sub>2</sub>O<sub>2</sub> /0.5% KOH medium until pigmentation has completely disappeared

*(Per ml: 100 ul of 30% solution of H<sub>2</sub>O<sub>2</sub> [this is how it comes], 250 ul of 2% KOH, 500 ul H<sub>2</sub>O).*

This takes between 30 min and 1 h; progress can be checked by observing the embryos under a dissecting scope. There will be a lot of bubbles formed in the medium – DO NOT CAP TUBES DURING THIS PROCESS – THEY WILL POP. (HYDROGEN PEROXIDE SOLUTION SHOULD BE PREPARED FRESH IMMEDIATELY BEFORE USE).

- b. Wash the embryos 2x for 5 min each in 1x PBS to remove the H<sub>2</sub>O<sub>2</sub> and stop the bleaching reaction.
  - c. Progressively dehydrate the embryos by washing for 5 min in 33% MeOH, 66% MeOH and 100% MeOH.
5. **Store in MEOH at -20°C for at least 2 hours before proceeding (or store long-term for later use).**

### Hybridization

How you start, depends on whether you are using freshly fixed embryos or embryos that have been stored in MeOH.

- Freshly fixed embryos:
  - a. Wash 5x, 5 min in PBT to remove PFA.
  - b. Remove chorions if not already done:
    - i. Transfer to Petri dishes and remove chorions.
    - ii. Transfer embryos back to microfuge tube (~10/tube).
    - iii. Wash 2x, 5 min PBT.
    - iv. Remove chorions, if not already done.
- MeOH embryos (from storage) must be rehydrated:
  - a. Wash 1X, 5 min in 66% MeOH:33% PBT.
  - b. Wash 1X, 5 min in 33% MeOH:66% PBT.
  - c. Wash 2X, 5 min each in PBT.

1. Treat with Proteinase K (ProtK) to permeabilize using the following incubation times at RT (or only if embryos are >20 somites).
  - a. Treat embryos with 10 µg/ml (in PBT) of Proteinase K for the following time periods. **DON'T OVER TREAT EMBRYOS!**
    - i. 3-5 somites            1 min
    - ii. 8-10 somites        3 min
    - iii. 16-20 somites      5 min
    - iv. 24 hpf                10 min
    - v. 30-36 hpf            15 min
  - b. Fix embryos again in 4% PFA for 20 min at RT.
  - c. Wash 5x, 5 min each with PBT.
2. Warm hybridization solution to 70°C.
3. Add ~500 µl warm prehyb solution to each tube with embryos.
- 4. Pre-hybridize embryos for 2-4 hours at 70°C.**
5. Prepare hybridization solution with labeled probe of choice – start with a dilution of 1:200 probe:hybridization solution in a volume of 200 µl per tube. New probes may need to be titered.
- 6. Hybridize O/N at 70°C in water bath.**

## Day 2

### Washing out the probe

1. **Prewarm** the first 5 wash solutions – washes.
2. Remove tubes from water bath and  
**SAVE THE PROBE IN THE HYBRIDIZATION SOLUTION AND STORE AT -20°C– IT CAN BE RE-USED!!**
3. 1x quick wash in prewarmed hyb solution
4. 1x, 5 min    66% hyb:33% 2X SSC at 65°C.
5. 1x, 5 min    33% hyb:66% 2X SSC at 65°C.
6. 1x, 5 min    2X SSC Tween at 65°C.
7. 2x, 20 min   0.2X SSC Tween at 65°C.
8. 1x, 5 min    66% 0.2X SSC:33% PBT at RT
9. 1x, 5 min    33% 0.2X SSC:66% PBT at RT
10. 1x, 5 min   PBT at RT

### Antibody (Ab) binding

1. Remove PBT and add 500 µl PI buffer. Incubate embryos for 2-4 hours in PI (blocking) buffer at RT to saturate non-specific binding sites for the antibody.
2. Prepare 1:5000 dilution of anti-dig Ab solution (1 µl anti-dig-AP Fab Ab to 5 ml of PI buffer). This dilution can be made using pre-bound Ab to further reduce non-specific staining.
  - a. Pre-bind Ab by making a 1:100 dilution of the antibody in PI buffer. Incubate with fixed/dechorionated embryos that have been rocked for at least 2 hours at RT. Pre-bound Ab can be stored at 4 °C. To make a final dilution of 1:5000 of the antibody – dilute 4 µl of pre-bound Ab per 200 µl of PI buffer)



3. Remove PI buffer and add ~200  $\mu$ l of 1:5000 dilution of anti-dig-AP Fab Ab.
4. Incubate **O/N at 4°C** with tubes on their sides (rocking is not necessary but make sure embryos are covered with Ab-containing buffer).

---

## Day 3

### Develop

1. Wash 2x rapidly with PBT.
2. Wash 6x, 15 min with PBT on rocker. (At this point you can store O/N at 4°C and wash the next day).
3. Wash 3x, 5 min in AP buffer. During washes, prepare developing solution (500  $\mu$ l per tube). Do not leave embryos in AP buffer for extended periods of time.
4. Remove wash solution and add 500  $\mu$ l of developing solution to each tube.
5. Transfer embryos to 9-well glass ISH dishes **USING A GLASS PASTEUR PIPET, (DO NOT USE A REGULAR PIPET TIP, THE EMBRYOS WILL STICK TO THE INSIDE WALL OF THE TIP!)**.
6. Store in the dark (place the top of a cardboard freezer box over the ISH dishes).
7. Monitor samples for blue product every 5-15 minutes.
8. Once desired intensity is reached, transfer embryos to clean microfuge tube.
9. Wash several times in stop solution.
10. Store in stop solution (in brown tubes) 4°C in the dark.

### Imaging embryos

1. Dehydrate embryos through a graded series of methanol solutions (30%, 50%, 70% and 100%) for 5 minutes each.
2. To clear the embryos, transfer them to the well of a spot plate containing a 1:2 solution of benzyl alcohol:benzyl benzoate (300-400  $\mu$ l per well). Allow embryos to sink in the well for 2-3 minutes.
3. Pick up the embryos with a glass hair (or the nylon loops) – NOT WITH A PIPET – try to reduce the volume of liquid transferred with the embryos! Place each embryo in a drop of Permount on a bridged coverslip\*.
4. Gently lower another coverslip on top of the embryo. Roll the upper coverslip until the embryo is in the desired orientation.

*\*Bridged coverslips can be made by gluing (with crazy glue), three ~22x22 coverslips on either end of a microscope slide.*




Figure 1. Embryos on left - without clearing. Embryos on right - with clearing.