

CELL CULTURE – GENERAL PROTOCOL

Media

225 ml DMEM (with 1.5 g/L sodium bicarbonate)
25 ml calf serum
2.5 ml Pen/strep
2.5 ml Glutamax
2.5 ml NEAA
2.5 ml sodium pyruvate
Filter

From frozen culture

1. Dilute vial of cells 1:10 in DMEM (with 10% calf serum)
2. Spin 5' @ setting 3 (Barnes lab)
3. Aspirate off media
4. Resuspend cells in 12 ml media
5. Plate on 100 mm cell culture dish

Splitting cells (1:4)

1. Aspirate off media
2. Rinse with 6 ml 1x filter-sterilized PBS and aspirate off PBS
3. Add 1 ml of trypsin (TripleE; Invitrogen) per 100 mm plate
4. Incubate at 37C for 3-5'
5. Inactivate trypsin with 9 ml complete media (i.e., w/serum) per 1 ml trypsin
6. Pipet up and down
7. Transfer to 15 ml conical bottom tube
8. Pipet up and down 15-20x in conical bottom tube
9. Spin down 2' at ~2000g
10. Aspirate off media – DON'T ASPIRATE THE PELLETT!
11. Add 10 ml media to 4 new plates
12. Resuspend cells in 8 ml
13. Pipet 2 ml per plate (for a total volume of 12 ml)

Freezing cells

1. Aspirate off media
2. Rinse with 6 ml 1x filter-sterilized PBS
3. Aspirate off PBS
4. Add 1 ml trypsin
5. Incubate at 37C for 3-5'
6. Inactivate trypsin with 9 ml of complete serum and transfer to a 15 ml conical bottom tube
7. Pipet up and down 15-20x in conical bottom tube
8. Spin down 2' at ~2000g
9. Aspirate off media
10. Add 1 ml serum/10%DMSO per 100 mm plate that was spun down
11. Pipet up and down and transfer to cryovials
12. Place vials in Styrofoam box in -80C
13. Transfer to liquid nitrogen tank after 24h at -80C

Transfection Protocol (12 well plate):

1. Day before transfection: plate cells at $0.4\text{--}2.0 \times 10^5$ cells/well (40-80% confluent)
 - a. **Make sure you have filtered medium that is free of proteins, antibiotics and serum for complex formation with Superfect reagent.**
2. On the day of transfection:
 - a. Dilute supercoiled DNA in TE buffer to 75 ul volume with cell growth medium containing no serum, proteins or antibiotics.
 - b. Spin down.
 - c. Add appropriate volume of Superfect reagent [7.5 ul] (Qiagen).
 - d. Incubate samples for 5-10 minutes at RT to allow complex formation.
 - e. While complex formation takes place, aspirate medium from dish and wash cells once with PBS.
 - f. Add appropriate volume of cell growth medium (containing serum and antibiotics – 400 ul) to the reaction tube containing transfection complexes. Mix by pipetting up and down 2x and transfer total volume to the cells.
 - g. Incubate cells with the transfection complexes for 2-3 hours under normal growth conditions
 - h. Remove medium by aspiration.
 - i. Wash cells 3-4x with PBS.
 - j. Add fresh complete cell growth medium.

Culture format	DNA:Superfect ratio	DNA (ug)	Final vol. of DNA in SF med.	Vol. Superfect reagent (ul)	Vol. of Serum-containing med (ul)
6-well	1:2	2	100	4	600
	1:5	2	100	10	600
	1:10	2	100	20	600
12-well	1:2	1.5	75	3	400
	1:5	1.5	75	7.5	400
	1:10	1.5	75	15	400

Use 12-well plate @ 1:5 ratio of DNA:Superfect to start and 1:1, 1:10 and 1:50 ratios of control:experimental constructs (control construct is phrl-CMV luciferin [from Promega and given to us by Daniel in Boyer's lab]):

1. pGL3-CMV (3 replicates each DMSO or TCDD)
2. wtFoxQ1b(-7.5)_Luc (3 replicates each DMSO or TCDD)
3. FoxQ1b(-7.5) Δ 1-7ahre_Luc (3 replicates each DMSO or TCDD)
4. 2dluc (from T. Gasiewicz @ U. Rochester) (3 replicates each DMSO or TCDD) – this construct has 2 consensus AhREs driving luciferase expression

Too many samples to start?

- alternatively #1-4 (3 replicates each DMSO or TCDD) at 1:10 ratios of control:experimental constructs.
- Total: 24 wells.

Plate reader set-up (for me)

1. KC4
2. System/reader/dispenser
 - a. Volume: 2 ml (70% EtOH in 50 ml conical tube); Rate: 300
 - b. Volume: 2 ml (H₂O); Rate 300
 - c. Prime air
 - d. Substrate (protected from light), prime #1 with firefly luciferase substrate (Volume: 1100) (make sure yellow is dispensed), prime #2 with Renilla luciferase substrate
 - e. Notes:
 - i. Need 5 ml for 4, 24-well plates
 - ii. Lyse in 100 ul (can use 50 ul)
 1. Rock 15' (or freeze O/N)
 2. Transfer 20 ul to 96 well plate
 3. (cleaner version – transfer to tube, vortex, spin down at ma for 5' and take 20 ul of supernatant)
 - f. Protocol:
 - i. Open dual luciferase
 - ii. Settings/layout/type sample
 - iii. Options/export/Microsoft excel
3. Clean
 - a. Purge reagent (1100 ul)
 - b. 70% EtOH
 - c. H₂O
 - d. air